



Guidelines for Ordering a SeaBrook Custom Abutment for Final Restoration

PLEASE FORWARD TO LABORATORY:

1. Full mouth impression + opposing (Non-alginate only, please)*
2. Bite relation (VPS registration material only)**
3. Copy of surgeon's post-implant info letter (see checklist below)
4. Lab slip with request for SeaBrook Custom Abutment & crown shade

NOTE: PLEASE ALLOW LAB A MINIMUM OF 2 WKS. TURN-AROUND FOR COMPLETION OF CASE. 3 WKS NEEDED IF GOLD-HUE COATING IS REQUESTED.

* Full mouth impression necessary because precision abutment design software factors in patient's full mouth anatomy. (Opposing model accepted in alginate but not preferred)

** VPS or Polyether material (vs. wax or alginate) provides the required accuracy needed for articulation and design of precision-machined abutments.

If thin tissue or all-ceramic crown case, please request gold-hue SeaBrook abutment.

OPTIONAL: Any additional pre-surgical models, surgical guides or cast of prosthetic device used for the case. These will all be optically scanned to assist in the design of your ideal custom abutment!

IMPLANT DATA CHECKLIST

1. Copy of Surgeon's Implant Info Sheet or please include on lab slip:

- Patient Name (lifetime warranty on abutment & implant for patient with SeaBrook abutment)
- Surgeon Name
- Tooth Number
- Implant Type
- Platform Diameter
- Healing Abutment/Collar Information (optional)

2. Please request an SeaBrook Abutment on the lab slip and what emergence profile width you are comfortable placing, based on the size of the healing collar (width preference is optional! If given no indication, SeaBrook will choose appropriate design preferences as in traditional method):

- **Large (Ideal Emergence)** – A releasing incision or gingivoplasty may be needed
- **Medium** – Blanching of soft tissue may be needed
- **Narrow** – Usually no/slight blanching and no incision required to place

Please call Tim Holbrook at (800) 693-7487 with questions.